

# capnography mounting options

for use with  **VIAMED** VM-2500 series of Capnographs



- 4410500** VM-2500-M - Central European
- 4410501** VM-2500-M - Asia
- 4410530** VM-2500-M - Veterinary



- 4410520** VM-2500-S - Central European
- 4410521** VM-2500-S - Asia
- 4410540** VM-2500-S - Veterinary



- 4410550** VM-2500-MG - Central European
- 4410560** VM-2500-MG - Veterinary



Universal Mounting Plate  
for VM-2500 Range  
**4420600** Rev. 1  
**4420601** Rev. 2



**0121180** V Wall Mount



**0121181** V Mount Pole Clamp - Small  
*Clamping Width: 14 -25mm*



**0121182** V Mount Pole Clamp - Large  
*Clamping Width: 16 - 40mm*



**0121184** V Mount Rail Clamp  
*Attaches to most types of  
30mm (H) x 10mm (D) Medirail*



**0121200** Universal Mounting Clamp  
*Attaches to round or flat objects up  
to 3.8 cm thick. Facilitates positioning  
of the device at any angle*

